





WECK HEM-O-LOK LIGATION SYSTEM

Safety without compromise





THE CLIP BANK

EASY LOADING WITH CONVENIENT ADHESIVE STRIP

- Our Hem-o-lok ligation system combines speed and safety in one single procedure.
- Choose from a variety of ligating clip sizes for the ligation of vessels and/or tissues between 2 and 16 mm.
- The Hem-o-lok clip bank allows for effortless and precise loading, thus making surgical interventions more efficient.
- Closure of the clip is confirmed by tactile and audible feedback.
- The colour code of the clip bank corresponds to that of the applier.
- Each cartridge contains 6 clips.
- A special adhesive strip at the bottom of each cartridge enables you to choose the best possible site for clip delivery.
- Unlike absorbable clips, Hem-o-lok ligating clips can be removed with our specially developed removers if necessary.
- The clips are inert, non-toxic and radiolucent.
- The cartridges adhere to any sterile surface.



THE APPLIERS

PERFORMANCE FEATURES AND CLEANING OPTIONS

- Thanks to the special clip design, closure of the clip can be clearly felt by the surgeon during application.
- Both curved and angled appliers are available for open surgery.
- Endoscopic Hem-o-lok appliers are available in sizes ML, L, and XL. They can be taken apart for easy cleaning, thus eliminating potential contamination.
- A 360° rotation knob facilitates access and improves visibility.
- The appliers are colour-coded to match the colours of the corresponding cartridges.

THE HEM-O-LOK LIGATION SYSTEM – PRECISION IN EVERY DETAIL

Teleflex is a leading manufacturer of metal and polymer clips. Our Weck Hem-o-lok ligation system forms a line of instruments perfectly adapted to each other, meeting all the special needs of surgeons, nurses and operating theatre personnel alike.



THE CLIP

TOP-CLASS FUNCTIONALITY

CLIP SIZE	LIGATION OF VESSELS AND TISSUES
● medium	2 – 7 mm
● medium-large	3 – 10 mm
● large	5 – 13 mm
● extra-large	7 – 16 mm

- Precision bosses engineered for maximum security.
- Penetrating lock: 1 x 2 teeth to provide effective closure.
- Bowed design enables added tissue to be encompassed.
- Tactile and audible response confirms that the Hem-o-lok ligating clip has locked onto vessel.
- Integrated teeth help prevent slippage in any direction.
- The surface properties of the non-absorbable polymer prevent infection.
- Living hinge provides flexibility when opening and closing.

LIGATING CLIPS



OPEN SURGERY



All Hem-o-lok ligating clips are sterile packaged. Ligating clip systems differ in their clip-closure characteristics. It is the responsibility of the user to confirm security of clips on vessels or other tissue. Not intended for contraceptive tubal occlusion.

HEM-O-LOK LIGATING CLIPS

WECK

REF.	CLIP SIZE	LIGATION OF VESSELS AND TISSUE	CLIPS/CARTRIDGE	CARTRIDGES/BOX	CLIPS/BOX
WK544220	● M	2 – 7 mm	6	14	84
WK544230	● ML	3 – 10 mm	6	14	84
WK544240	● L	5 – 13 mm	6	14	84
WK544250	● XL	7 – 16 mm	6	14	84

APPLIERS

WECK

REF.	CLIP SIZE	WORKING LENGTH	QTY
curved			
WK544113	● M	20 cm	1
WK544170	● ML	20 cm	1
WK544180	● L	20 cm	1

APPLIERS, LONG

WECK

REF.	CLIP SIZE	WORKING LENGTH	QTY
curved			
WK544115	● M	27 cm	1
WK544171	● ML	27 cm	1
WK544181	● L	27 cm	1
WK544191	● XL	27 cm	1

APPLIERS, LONG, 70° ANGLED

WECK

REF.	CLIP SIZE	WORKING LENGTH	QTY
WK544114	● M	27 cm	1
WK544172	● ML	27 cm	1
WK544179	● L	27 cm	1
WK544192	● XL	27 cm	1

CLIP REMOVERS

WECK

REF.	CLIP SIZE	WORKING LENGTH	QTY
WK544123	● M	27 cm	1
	● ML		
	● L		
WK544124	● L	27 cm	1
	● XL		

● M medium ● ML medium-large ● L large ● XL extra-large

ENDOSCOPIC SURGERY



HEM-O-LOK TAKE-APART APPLIERS

WECK

REF.	CLIP SIZE	Ø APPLIER	Ø TROCAR	WORKING LENGTH	QTY
straight					
WK544965T	● ML	5 mm	5 mm	33 cm	1
curved					
WK544995T	● L	10 mm	10 mm	33 cm	1
WK544990T	● XL	10 mm	10 or 12 mm	33 cm	1

HEM-O-LOK TAKE-APART APPLIERS, 20° ANGLED

WECK

REF.	CLIP SIZE	Ø APPLIER	Ø TROCAR	Ø SILS	WORKING LENGTH	QTY
WK544965T20	● ML	5 mm	10 mm	5 mm	33 cm	1
WK544995T20	● L	10 mm	12 mm	10 mm	33 cm	1
WK544990T20	● XL	10 mm	15 mm	12 mm	33 cm	1

HEM-O-LOK NON-TAKE-APART APPLIER, LONG

WECK

REF.	CLIP SIZE	Ø APPLIER	Ø TROCAR	WORKING LENGTH	QTY
straight					
WK544965L	● ML	5 mm	5 mm	45 cm	1
curved					
WK544995L	● L	10 mm	10 mm	45 cm	1
WK544990L	● XL	10 mm	10 or 12 mm	45 cm	1

HEM-O-LOK AUTO ENDO 5 APPLIERS

WECK

REF.	CLIP SIZE	Ø APPLIER	Ø TROCAR	WORKING LENGTH	QTY
straight					
WK543965	● ML	5 mm	5 mm	33 cm	3

CLIP REMOVERS

WECK

REF.	CLIP SIZE	Ø REMOVER	Ø TROCAR	Ø SILS	WORKING LENGTH	QTY
straight						
WK544121T	● ML	5 mm	5 mm	–	33 cm	1
WK544121L	● ML	5 mm	5 mm	–	45 cm	1
curved						
WK544130T	● L/XL	10 mm	10 or 12 mm	–	33 cm	1
WK544130L	● L/XL	10 mm	10 or 12 mm	–	45 cm	1
20° angled						
WK544130T20	● L/XL	10 mm	12 or 15 mm	12 mm	33 cm	1

ENDOSCOPIC TAKE-APART APPLIER

Our take-apart endoscopic appliers can be detached for easy cleaning. The 360° knob facilitates access and improves visibility.



Teleflex is a leading global provider of specialty medical devices used for diagnostic and therapeutic procedures in critical care, urology and surgery. Our mission is to provide solutions that enable healthcare providers to improve outcomes and enhance patient and provider safety. We specialise in devices for general and regional anaesthesia, cardiac care, respiratory care, urology, vascular access and surgery and we serve healthcare providers in more than 150 countries. Teleflex also provides specialty products for medical device manufacturers.

Our well known brands include ARROW®, DEKNATEL®, GIBECK®, HUDSON RCI®, KMEDIC®, LMA™, PILLING®, PLEUR-EVAC®, RÜSCH®, SHERIDAN®, TAUT®, TFX OEM®, VASONOVA®, VIDACARE® and WECK®, all of which are trademarks or registered trademarks (in the U.S. and/or other countries) of Teleflex Incorporated.

Teleflex global operations: Australia, Austria, Belgium, Canada, China, Czech Republic, France, Germany, Greece, India, Ireland, Italy, Japan, Korea, Malaysia, Mexico, Netherlands, New Zealand, Portugal, Singapore, Slovak Republic, South Africa, Spain, Switzerland, United Kingdom, Uruguay and USA.

YOUR INTERNATIONAL CONTACTS:

TELEFLEX HEADQUARTERS INTERNATIONAL, IRELAND

Teleflex Medical Europe Ltd., IDA Business and Technology Park,
Dublin Road, Athlone, Co Westmeath
Phone +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73
orders.intl@teleflex.com

AUSTRALIA/NEW ZEALAND 1300 360 226

AUSTRIA +43 (0)1 402 47 72

BELGIUM +32 (0)2 333 24 60

CHINA (SHANGHAI) +86 (0)21 6163 0965

CHINA (BEIJING) +86 (0)10 6418 5699

CZECH REPUBLIC +420 (0)495 759 111

FRANCE +33 (0)5 62 18 79 40

GERMANY +49 (0)7151 406 0

GREECE +30 210 67 77 717

INDIA +91 (0)44-2836 5040

ITALY +39 0362 58 911

JAPAN +81 (0)3 3379 1511

KOREA +82 2 536 7550

NETHERLANDS +31 (0)88 00 215 00

PORTUGAL +351 22 541 90 85

SINGAPORE +65 6439 3000

SLOVAK REPUBLIC +421 (0)3377 254 28

SOUTH AFRICA +27 (0)11 807 4887

SPAIN +34 918 300 451

SWITZERLAND +41 (0)31 818 40 90

UNITED KINGDOM +44 (0)1494 53 27 61

For detailed information see www.teleflex.com

Teleflex is a registered trademark of Teleflex Incorporated or its affiliates.
The products in this catalogue may not be available in all countries. Please contact your local representative.
All data current at time of printing (05/2016). Subject to technical changes without further notice.
© 2016 Teleflex Incorporated. All rights reserved.

94 05 70 - 00 00 01 · REV A · MC / WM · 05 16 01