





XPS Nexus™
ENT Power Console

**TAKING
SINUS SURGERY
FURTHER, TOGETHER**



Medtronic

XPS NEXUS™ ENT SYSTEM



KEY FEATURES



Connects the M2, M4, and Indigo™



Handpiece recognition



Portable, easy to transport, and lightweight



IV pole (*sold separately*)



Built-in irrigation pump



Intuitive design with a user-friendly interface



Quick and simple set up



Foot pedal with non-slip design (*sold separately*)

A Pump Panel

1. Flow selection (CC/min)
2. Flow increase

B M2/M4 Panel

3. Oscillation mode
4. RPM selection
5. Forward mode
6. RPM increase/decrease

C Indigo Panel

7. Forward
8. Reverse
9. RPM selection
10. RPM increase/decrease

D Irrigation Pump

E Power/Connection Panel

11. Indigo, M4/M2 receptacle
12. Power light (green)
13. Error indication light (amber)
14. Handpiece alignment mark
15. Foot pedal receptacle

HANDPIECES AND ATTACHMENTS

Handpiece Features

StraightShot™ M4 Microdebrider

- Rotating tip curved blade design allows 360-degree cutting operation
- Ergonomic handpiece with rotating finger wheel facilitates simple, intraoperative blade tip rotation
- 60-degree angled blade provides access to frontal recess

StraightShot™ Magnum II Microdebrider

- Operates at speeds up to 12,000 RPM (*forward*) and 5,000 RPM (*oscillating*)

Indigo™ High-Speed Otologic Drill

- High speed and ample torque for fast debulking of the mastoid
- Available with straight and angled attachments
- Improved surgical access and visibility in narrow otologic spaces
- Works with wide variety of Medtronic Xomed otology burs
- Grooves on attachment nose design allow for more stable irrigation

StraightShot™ M4



StraightShot™ Magnum™ II



Indigo™ High-Speed Otologic Drill



ORDERING INFORMATION

Description	Part Number
XPS Nexus™ Console.....	1894510
XPS Nexus™ Foot Control	1894520
XPS Nexus™ IV Pole	1894515
XPS Nexus™ Carrying Case.....	1894525
Straightshot™ M4 Microdebrider	1898200T
StraightShot™ M4 Instrument Tray	1898075
StraightShot™ Magnum II Microdebrider	1897200
StraightShot™ Magnum II Instrument Tray	1896840
Indigo™ High-Speed Otologic Drill	1845000
Indigo™ Otologic Straight Attachment	1845010
Indigo™ Otologic Angled Attachment	1845020
Indigo™ Instrument Tray.....	1845040
Indigo™ Tubing	1845030
IPC™ System Power Cord , Standard, North America, 3M.....	1895820
IPC™ System Power Cord , Standard, North America, 6M.....	1895821
IPC™ System Power Cord , UK/Ireland, 6M.....	1895824
IPC™ System Power Cord , Continental Europe, 6M.....	1895825
IPC™ System Power Cord , Japan	1895823



Rx only. Refer to product instruction manual/package insert for instructions, warnings, precautions and contraindications.

For further information, please call Medtronic ENT at 800.874.5797 or consult Medtronic's website at www.medtronic.com/ent.

Medtronic

Medtronic plc ENT

6743 Southpoint Drive N
Jacksonville, FL 32216
USA
Toll free: (800) 874-5797
Telephone: (904) 296-9600
Fax: (800) 678-3995

International Telephone Numbers

Adriatic Region 385-1-488-1120	Israel 972-9-972-4400
Argentina 54-11-4898-5700	Italy 39-02-24137-324
Australia 1-800-668-670	Japan 81-3-6430-2017
Baltic Region 37-1-67560226	Korea 82-2-3404-3600
Belgium 32-2456-09-09	Lebanon 961-1-370-670
Brazil 55-11-2182-9200	Luxembourg 32-2456-09-09
Canada 1-800-217-1617	Malaysia 60-37-953-4800
Chile 56-2-2655-5110	Mexico 52-55-11-02-90-30
China 86-21-20325888	Netherlands 31-45-566-8800
Colombia 57-1-742-7300	Poland 48-22-4656900
Czech Republic 420-2-9657-9580	Russian Fed. 7-495-580-73-77
France 33-470-679-800	Singapore 65-6776-6255
Germany 49-2159-8149-209	South Africa 27-11-466-1820
Greece 30-210-67-79-099	Spain 34-91-625-05-40
Hong Kong 852-2919-1300	Taiwan 886-2-2183-6000
Hungary 36-188-90600	Thailand 662-232-7400
India 91-22-26836733	UK 44-1923-205-166